

GLISSANDO IMAGESHARE

WHOLE SLIDE IMAGE ARCHIVING AND SHARING SOLUTION



EASILY ADD WSI STORAGE AND EFFICIENT REMOTE VIEWING TO YOUR GLISSANDO SCANNER

The Glissando ImageShare is a complete, stand-alone solution for storing and sharing large whole slide image (WSI) scans. The ImageShare workstation links directly to one or more Glissando slide scanners, providing virtually immediate availability for WSI viewing both locally and remotely.

The ImageShare's innovative direct link sharing software allows users to easily share WSI results via efficient browser-based viewing of folders and individual WSI results on the ImageShare, removing the time-consuming step of uploading large WSI images to a secondary server or cloud solution.

When combined with the Glissando medium- and low-capacity slide scanners, the Glissando ImageShare provides a perfect solution for quickly and easily implementing whole slide image scanning, storage, and remote sharing into your everyday workflows.

KEY FEATURES

- ◆ **Easy-to-implement, highly efficient remote WSI viewing**
- ◆ **Direct connection to one or more Glissando scanners**
- ◆ **Stand-alone solution needs only an internet connection**
- ◆ **No external upload of large WSI scans required**
- ◆ **60 terabytes of archiving space**
- ◆ **Shareable browser links directly to folders and individual images**
- ◆ **Fault-tolerant, hot swappable RAID disk storage array**

IMAGESHARE SOLUTION OVERVIEW

GLISSANDO
SCANNER(S)

GLISSANDO
IMAGE SHARE

IMAGESHARE WSI
FOLDER/IMAGE

INTERNET LINK TO
FOLDER/IMAGE

REMOTE VIEWING
VIA BROWSER



SPECIFICATIONS

Scanner connection	Direct Ethernet link between scanner and ImageShare
Remote connection	High-speed internet
Remote viewer architecture	Browser-based
Remote view options	Direct browser link to a specific folder or image
WSI storage capacity	60 TB
Supported WSI format	SVS
System OS	Microsoft Windows® 10 IoT Enterprise LTSCB
System Ethernet ports	2
System wireless	Yes, Wi-Fi included, 4G/5G optional*

*4G/5G option requires 4G/5G provider plan not-included

THE OBJECTIVE IMAGING DIGITAL PATHOLOGY FAMILY

GLISSANDO



Compact, high quality scanner with a capacity for 2 slides.

GLISSANDO 20SL



Auto-loading 20-slide scanner for medium capacity scenarios.

GLISSANDO POL



Scanning in both brightfield and polarised light.

GLISSANDO
IMAGE SHARE



Stand-alone whole slide image archiving and sharing workstation.



Objective Imaging Ltd.
Unit 1, Quay Court, Colliers Lane, Stow cum Quy
Cambridge, CB25 9AU
United Kingdom



+44 1223 813777



sales@objectiveimaging.com



www.objectiveimaging.com



FM631706